



EST 2002

lipoelastic[®]

MEDICAL PRODUCTS

lipoelastic book

POSTOPERATIVE COMPRESSION GARMENTS

CONTENT

THE STORY OF LIPOELASTIC	5
BENEFITS	6
COMPRESSION THERAPY	10
HOW DO I MEASURE MYSELF?	13
HOW WE DEFINE V AND T GARMENTS?	14
TYPES OF FASTENING	17
TYPES OF FINISHING	17
TYPES OF HYGIENE OPENING FOR WOMEN	18
BREAST GARMENTS AND STABILIZER BANDS	20
PI elite	21
PI special	22
PI super	23
PI perfect	24
PI plus	25
PI ideal	26
PI active	27
PI relax	28
PI extra	29
PI filling	30
LIPOBRA	31
PS ideal	32
PS special	33
SG	34
SI 01	35
SI formed	36
PI unique	37
PI premium	38
FEMALE GARMENTS	40
VH Comfort	41
VH	42

VH Variant	43
VH low Variant	44
VH unique Variant	45
VH body Variant	46
VH special Comfort	47
VH special	48
VF Comfort	49
VF Variant	50
VF body Variant	51
VF body BBO Variant	52
VF body BBC Variant	53
VD Comfort	54
VD Variant	55
VD body Variant	56
VD body BBO Variant	57
VD body BBC Variant	58
VD special Comfort	59
VB Comfort	60
VB Variant	61
VBfLg Variant	62
TF Comfort	63
TF Variant	64
TF leggings	65
TD Comfort	66
TD Variant	67
TD leggings	68
TB Comfort	69
TB Variant	70
TB leggings	71
TBfLg Variant	72

TBfLg leggings	73
ACTIVE leggings	74
ARM SLEEVES	76
AP Variant	77
AP long Variant	78
APf Variant	79
AS Variant	80
AS long Variant	81
FEMALE VESTS AND BODY SUITS	82
MT smooth Variant	83
MTf smooth Comfort	84
MH Comfort	85
MH special Comfort	86
MHF Comfort	87
MHD Comfort	88
MGF Variant	89
MGD Variant	90
MALE GARMENTS	92
VHmS Comfort	93
VHmM Comfort	94
VFm Comfort	95
TDm Variant	96
MTmS Comfort	97
MTm Comfort	98
MTmL Comfort	99
MTmL Variant	100
MGmm Comfort	101
MGm long Variant	102
ABDOMINAL BINDERS AND ACCESSORIES	104
KPress	105

KPlus	106
KP soft	107
KP special	108
KPG	109
KPad	110
KPad 01	110
KPad 02	111
KPad 03	111
KPad 04	112
KPad 4D	112
FACIAL GARMENTS	113
FM 02	114
FM 03	115
FM special	116
FM extra	117
FM 01	118
PU 03	119
OTHER GARMENTS	120
LIPO slide	121
LIPO extend T and V	122
LIPO extend AP	122
ANTI-EMBOLISM STOCKINGS	123
LIPOTHROMBO AD	124
LIPOTHROMBO AG	125
SCAR TREATMENT	126
LIPOELASTIC sheet LOLLIPOP	127
LIPOELASTIC sheet ANCHOR	127
LIPOELASTIC sheet STRIP01	128
LIPOELASTIC sheet STRIP02	128
LIPOELASTIC gel	129



“We have always believed that employees are, in fact, the company’s shareholders: they form its foundation and influence its destiny. Therefore, we consider them part of our LIPOELASTIC family; a family with a global reach.”

Alena and Tomáš Hradil,
founders of LIPOELASTIC a.s.

THE STORY OF LIPOELASTIC

We are LIPOELASTIC a.s., a Czech family company founded in 2002 by Alena and Tomáš Hradil.

From a small company with 15 employees, we have developed into one of the world's leading manufacturers of post-operative compression garments. Today, more than 20 years later, our company has over 220 employees, and we continue to grow. Thanks to the work ethic of our dynamic team, we are currently exporting to more than 80 countries, and we are among the leading providers of compression garments in the world.

Present since

2002

ON THE MARKET

Wide portfolio of

100

PRODUCTS

Export to

80

COUNTRIES

BENEFITS



R&D Department

Despite the success of our products, we continue to look for new ways to improve. Our R&D department works hard to constantly innovate and develop our products, and make the whole team's work even better. Thanks to our cooperation with the best surgeons worldwide, we keep improving our garments to enable a perfect fit and maximum comfort. The greatest reward for us is a satisfied patient.



Modern Technologies

We constantly innovate and use the most modern technological procedures in our production. For example, ultrasonic and hemming technologies enable the joining of textiles by welding and gluing, making our products highly durable and less prone to rip.



Our Certifications

Quality and safety are our top priorities. Our products meet all the qualitative and material standards incorporated in the following certifications:



OEKO-TEX® STANDARD 100

LIPOELASTIC® products are awarded the OEKO-TEX® STANDARD 100 certificate, 02.0.8039 Hohenstein HTTI, which guarantees the absence of harmful and is therefore harmless for human health.



ISO 13485 Certification

LIPOELASTIC® products are manufactured in a facility whose quality management system is certified by the ISO 13485 standard (ČSN EN ISO 13485: 2016).



Declaration of Conformity

LIPOELASTIC® products are in conformity with the current EU legal regulations concerning medical devices and hold the CE marking.



Transportation of Sweat

Optimal values of sweat transfer in fabrics are considered to be in the range of 800 – 900 g/m²h. The LIPOELASTIC® materials reach values of up to 815 g/m²h.



Moisture Removal

Tests performed by the Czech Textile Testing Institute have proven the moisture removal and breathability of LIPOELASTIC materials.

BENEFITS



The Precision of Our Products

We do not compromise when it comes to quality of our products. All LIPOELASTIC® garments are handmade by our team of experienced sewists, and each post-operative garment is inspected in detail before dispatch from our factory.



Exclusive Materials

Materials for our products are carefully selected to meet our quality standards. We exclusively buy the materials of the highest quality and only from our trusted suppliers in the European Union.



Flat Seams

The comfort of LIPOELASTIC® post-operative compression garments is determined by small details that can make all the difference. One of these would be our unique flat seams technology ensuring maximum comfort.



Seamless Moulded Cups

The breasts are very sensitive after plastic surgery. We shape the cups of most of our bras by heat to prevent seams in sensitive areas.



Cotton Bras

Cotton is a natural and breathable material. For the production of some of our post-operative bras, we have chosen a material with a high percentage of cotton to provide patients with the greatest comfort and best post-operative care.



Made-to-Measure

We understand that every patient is unique and has different needs. Therefore, at LIPOELASTIC, we also manufacture made-to-measure garments based on a request from the patient's surgeon.



Textile Label

Placing the label on the outside of the garment prevents unpleasant scratching. The textile labels are also colour-matched to the colour of the garment.



Less Packaging

Join our journey towards sustainability by implementing new storage solutions and reduced packaging options.


COMPRESSION THERAPY

Compression therapy is an essential part of the recovery process following most plastic and aesthetic procedures. The application of compression pressure promotes the circulation of bodily fluids, effectively reducing their accumulation in the surgical area.*

Wearing LIPOELASTIC® compression garments after plastic surgery enhances patients' comfort during their recovery. We place high standards on the design and testing of the input parameters of our materials to achieve the recommended compression.


**Dixit VV, Wagh MS. Unfavourable outcomes of liposuction and their management. Indian J Plast Surg. 2013;46(2):377-392. doi: 10.4103/0970-0358.118617*

**Skin without
compression**



Illustrative picture

**Skin with
compression garment**



Illustrative picture





THE HEALING POWER OF COMPRESSION

Minimises swelling
and bruising



Pain
management



Prevents
complications



Tissue
reshaping



Scar
management



Speeds up the
recovery process



HOW TO MEASURE CORRECTLY?

Choosing the right-sized postoperative garments is necessary to achieve the right compression level. Therefore, it is essential to pay extra attention while measuring the required body parts.

- 1** *For the best results, it is recommended to always measure naked. Stand upright with a relaxed belly, and breathe naturally.*
- 2** *The measuring tape must be in a horizontal position. It should not point downwards or upwards to ensure that the measurement is not distorted.*
- 3** *When taking the measurements, the measuring tape should not be too tight or too loose. It should come to contact with the body only lightly.*
- 4** *Write down the measured values, and then re-measure to double-check for accuracy.*
- 5** *Using your measurements, choose the correct size of the garment in line with the product size chart.*

Are You In Between Two Sizes?

If you find that your measured size is in between two sizes in the product charts, we always recommend that you choose the smaller size to ensure sufficient compression. In cases where the size cannot be determined exactly, we offer the option of custom-made products. Custom-made production is only possible with compression garments.



Our TIP: Garments with Variant fastening (with 3 rows of hooks and eyes) can be adjusted as any post-operative swelling reduces.

HOW DO I MEASURE MYSELF?



Please measure the chest circumference for male garments.

HOW WE DEFINE V AND T GARMENTS?

LIPOELASTIC® compression garments are named according to the letters corresponding with a particular body part.

V type garments

Starting point is under bust



T type garments

Starting point is waistline




From the indicated points **V** and **T**, we determine the starting point of the garment. The length of the garment is indicated by points **H** to **B**.

The diagram on the right highlights our **VB** garment. As you can see, point **V** marks the starting point of the garment, and point **B** marks the finish.



TIP: We recommend opting for a garment which is almost 10 - 15 cm below the treated area to avoid any possibilities of excessive swelling.



NAME EXAMPLES OF OUR GARMENTS WITH TROUSERS





TYPES OF FASTENING

- 1**  **COMFORT** Zip fastening supported with a row of hooks and eyes to ensure easy and hassle-free use immediately after the procedure.
- 2**  **VARIANT** Adjustable fastening with 3 rows of hooks and eyes enables regulation of size and degree of compression as post-operative swelling reduces.
- 3**  **WITHOUT zipper**

TYPES OF FINISHING

Hemming



The hemming technology enables joining of textiles by gluing. An elegant hemming finish fully replaces the elastic band while maintaining its compression and elasticity.

Ultrasonic



Ultrasonic technology enables joining of textiles by welding. Finishing with this technology is very durable and the fabric is less prone to tear.

Lace



Contemporary lace without a silicone tape does not cut into the legs, is sufficiently fixed in the area and does not roll up.

Clean-cut



Certain types of garments are finished only with a clean-cut preventing the garment from rolling up or pressing in any way.

TYPES OF HYGIENE OPENING FOR WOMEN



Open crotch

Female garments with pants have an open crotch for greater comfort and accessibility.



Adjustable hygiene opening

Female garments without pants have a hygiene opening with three adjustable rows of hooks and eyes.



Without hygiene opening

Leggings are made without a hygiene opening.

TYPES OF HYGIENE OPENING FOR MEN



Hygiene opening - crotch

Male panty garments without pants have a high cotton content in the crotch area for maximum comfort.



Open crotch

Male trouser garments have an open crotch for better accessibility.



For maximum comfort of sensitive areas after plastic surgery, some LIPOELASTIC® garments are reinforced with a double padding in the abdomen or thigh areas.



EST 2002

lipoelastic®


MEDICAL PRODUCTS

LIPOELASTIC

BREAST GARMENTS
AND STABILIZER
BANDS

PI elite

LIPO-PI33Z00X



Use: Postoperative care after plastic-aesthetic procedures: Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors: ●

Product description:

- 1** soft comfortable underlay with zip protection
- 2** zip with hook and eye closure underneath for easy closing
- 3** clean-cut comfort finish
- 4** easy to use adjustable straps
- 5** seamless molded cups

NEW




Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Breast circumference	
		Cup B-C	Cup D-E
65	63-67	79-83	83-87
70	68-72	84-88	88-92
75	73-77	89-93	93-97
80	78-82	94-98	98-102
85	83-87	99-103	103-107
90	88-92	104-108	108-112
95	93-97	109-113	118-122



PI special

LIPO-PI01C00X



Use: Postoperative care after plastic-aesthetic procedures: Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors:  



COTTON

Product description:

- 1** straps with 4 positions of hooks-and-eyes closure
- 2** front zipper underlain with hooks
- 3** moulded cups reinforced with a double layer of fabric for perfect comfort and support
- 4** soft elastic band with logo

BEST
SELLER




Size	MEASUREMENTS IN CM / EU SIZE					
	Under breasts	Breast circumference				
		Cup A	Cup B	Cup C	Cup D	Cup E
60	58-62	72-74	74-76	76-78	78-80	-
65	63-67	77-79	79-81	81-83	83-85	-
70	68-72	82-84	84-86	86-88	88-90	90-92
75	73-77	87-89	89-91	91-93	93-95	95-97
80	78-82	92-94	94-96	96-98	98-100	100-102
85	83-87	97-99	99-101	101-103	103-105	105-107
90	88-92	102-104	104-106	106-108	108-110	110 - 112
95	93-97	107-109	109-111	111-113	113-115	115 - 117
100	98-102	-	114-116	116-118	118-120	120 - 122



PI super

LIPO-PI99V00X



Use: Postoperative care after plastic-aesthetic procedures: Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors: ○ ●



COTTON

Product description:

- 1** straps with 4 positions of hooks-and-eyes closure
- 2** 3 positions of hooks-and-eyes closure
- 3** moulded cups reinforced with a double layer of fabric for perfect comfort and support
- 4** soft elastic band with logo




Size	MEASUREMENTS IN CM / EU SIZE	
	Under breasts	Breast circumference
XS	63-67	77-81
S	68-72	81-88
M	73-78	88-97
L	79-84	94-101
XL	85-90	101-108
XXL	91-96	107-116
3XL	97-110	113-122



PI perfect

LIPO-PI23V00R




Use: Postoperative care after plastic-aesthetic procedures: Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors: 

Product description:

- 1** straps with 4 positions of hooks-and-eyes closure
- 2** 3 positions of pink hooks-and-eyes closure
- 3** moulded cups
- 4** seamless band finished with hemming technology




Size	MEASUREMENTS IN CM / EU SIZE					
	Under breasts	Breast circumference				
		Cup A	Cup B	Cup C	Cup D	Cup E
65	63-67	77-79	79-81	81-83	83-85	-
70	68-72	82-84	84-86	86-88	88-90	90-92
75	73-77	87-89	89-91	91-93	93-95	95-97
80	78-82	92-94	94-96	96-98	98-100	100-102
85	83-87	97-99	99-101	101-103	103-105	105-107
90	88-92	102-104	104-106	106-108	108-110	-
95	93-97	107-109	109-111	111-113	113-115	-
100	98-102	-	114-116	116-118	118-120	-



PI plus

LIPO-PI19V00R



Use: Postoperative care after plastic-aesthetic procedures: Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors: 

Product description:

- 1** straps with 4 positions of hooks-and-eyes closure
- 2** 3 positions of hooks-and-eyes closure
- 3** moulded cups
- 4** double underlay of fabric in the back for increased support
- 5** seamless band finished with hemming technology




Size	MEASUREMENTS IN CM / EU SIZE	
	Under breasts	Breast circumference
XS	63-67	77-81
S	68-72	81-88
M	73-78	88-97
L	79-84	94-101
XL	85-90	101-108
XXL	91-96	107-116
3XL	97-110	113-122



PI ideal

LIPO-PI04V00R



Use: Postoperative care after plastic-aesthetic procedures: Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors: 

Product description:

- 1** straps with 4 positions of hooks-and-eyes closure
- 2** 3 positions of hooks-and-eyes closure
- 3** moulded cups
- 4** soft elastic band with logo




Size	MEASUREMENTS IN CM / EU SIZE	
	Under breasts	Breast circumference
XS	63-67	77-81
S	68-72	81-88
M	73-78	88-97
L	79-84	94-101
XL	85-90	101-108
XXL	91-96	107-116
3XL	97-110	113-122



PI active

LIPO-PI03V00X



Use: Postoperative care after plastic-aesthetic procedures: Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors: ○ ●



COTTON

Product description:

- 1** straps with an adjustable buckle
- 2** 3 positions of hooks-and-eyes closure
- 3** moulded cups
- 4** soft elastic band with logo




Size	MEASUREMENTS IN CM / EU SIZE	
	Under breasts	Breast circumference
XS	63-67	77-81
S	68-72	81-88
M	73-78	88-97
L	79-84	94-101
XL	85-90	101-108
XXL	91-96	107-116
3XL	97-110	113-122



PI relax

LIPO-PI26V00X



Use: Postoperative care after plastic-aesthetic procedures: Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors:  

Product description:

- 1** adjustable tight straps
- 2** 3 positions of hooks-and-eyes closure
- 3** seamless band finished with hemming technology
- 4** moulded cups




Size	MEASUREMENTS IN CM / EU SIZE	
	Under breasts	Breast circumference
XS	63-67	77-81
S	68-72	81-88
M	73-78	88-97
L	79-84	94-101
XL	85-90	101-108
XXL	91-96	107-116
3XL	97-110	113-122



PI extra

LIPO-PI02V00X



Use: Postoperative care after plastic-aesthetic procedures: Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors: ○ ●

Product description:

- 1** straps with 4 positions of hooks-and-eyes closure
- 2** 3 positions of hooks-and-eyes closure
- 3** soft elastic band with logo




Size	MEASUREMENTS IN CM / EU SIZE	
	Under breasts	Breast circumference
XS	63-67	77-81
S	68-72	81-88
M	73-78	88-97
L	79-84	94-101
XL	85-90	101-108
XXL	91-96	107-116
3XL	97-110	113-122



PI filling

LIPO-PI21V00X



Use: Postoperative care after breast augmentation using the autologous fat transfer method

Colors: ●

Product description:


- 1** 3 rows of hooks-and-eyes closure in the back
- 2** fine net
- 3** seamless band finished with hemming technology

Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Breast circumference	
		Cup A - B	Cup C - D
65	63-67	77-81	81-85
70	68-72	82-86	86-90
75	73-77	87-91	91-95
80	78-82	92-96	96-100
85	83-87	97-101	101-105
90	88-92	102-106	106-110
95	93-97	107-111	111-115



LIPOBRA

LIPOBRA®



Use: Postoperative care after plastic-aesthetic procedures: Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy, Mastectomy, Breast Ablation

Colors: ●

Product description:

- 1** two pockets for the prosthesis
- 2** seamless cups
- 3** the cups are reinforced with a double layer of fabric for optimal support
- 4** adjustable shoulder straps
- 5** front hook and eye opening


S M L XL XXL

Under breast in cm	Breast circumference					
68-72	80-82	82-84	84-86	86-88	88-90	90-92
73-77	85-87	87-89	89-91	91-93	93-95	95-97
78-82	90-92	92-94	94-96	96-98	98-100	100-102
83-87	95-97	97-99	99-101	101-103	103-105	105-107
88-92	100-102	102-104	104-106	106-108	108-110	110-112
93-97	105-107	107-109	109-111	111-113	113-115	115-117
98-102	-	112-114	114-116	116-118	118-120	120-122
103-107	-	117-119	119-121	121-123	123-125	125-127
108-112	-	122-124	124-126	126-128	128-130	-
113-117	-	127-129	129-131	-	-	-



PS ideal

LIPO-PS04V00X



Use: Postoperative care after breast augmentation

Implant stabilizer band with compression bra as a stabilization support of implants following augmentation.

Colors: 

Product description:

- 1** straps with 4 positions of hooks-and-eyes closure
- 2** non-detachable implant stabilizer band with velcro fastening
- 3** breathable structure
- 4** moulded cups
- 5** 3 positions of hooks-and-eyes closure
- 6** seamless band finished with hemming technology




Size	MEASUREMENTS IN CM / EU SIZE	
	Under breasts	Breast circumference
XS	63-67	77-81
S	68-72	81-88
M	73-78	88-97
L	79-84	94-101
XL	85-90	101-108
XXL	91-96	107-116
3XL	97-110	113-122



PS special

LIPO-PS01C00X



Use: Postoperative care after breast augmentation

Implant stabilizer band with compression bra as a stabilization support of implants following augmentation.


Colors: ○ ●



COTTON

Product description:

- 1** straps with 4 positions of hooks-and-eyes closure
- 2** PI special bra in set with implant stabilizer band SI formed
- 3** zipper underlain with hooks
- 4** moulded cups reinforced with a double layer of fabric for perfect comfort and support
- 5** soft elastic band with logo
- 6** removable ergonomically shaped breast stabilizer




Size	MEASUREMENTS IN CM / EU SIZE				
	Under breasts	Breast circumference			
		Cup A	Cup B	Cup C	Cup D
65	63-67	77-79	79-81	81-83	83-85
70	68-72	82-84	84-86	86-88	88-90
75	73-77	87-89	89-91	91-93	93-95
80	78-82	92-94	94-96	96-98	98-100
85	83-87	97-99	99-101	101-103	103-105
90	88-92	102-104	104-106	106-108	108-110



SG

LIPO-SG00F00X



Use: Postoperative care after breast augmentation

Implant stabilizer band for separate use or in combination with surgical bra as a stabilization support of implants following augmentation.

Colors: ○ ●

Product description:

- 1** ergonomically shaped stabilizer band
- 2** velcro fastening
- 3** stronger, dense material - developed for larger implants
- 4** breathable material




MEASUREMENTS IN CM / EU SIZE	
Size	Breast circumference
XS	< 86
S	86-92
M	92-98
L	98-108
XL	108-114
XXL	114-122



SI 01

LIPO-SI00F00X



Use: Postoperative care after breast augmentation

Implant stabilizer band for separate use or in combination with surgical bra as a stabilization support of implants following augmentation.

Colors: ○ ●

Product description:

- 1** straight cut
- 2** velcro fastening
- 3** breathable material




	MEASUREMENTS IN CM / EU SIZE
Size	Breast circumference
XS/S	72 - 87
M/L	88 - 103
XL/XXL	104 - 120



SI formed

LIPO-SI09F00X



Use: Postoperative care after breast augmentation

Implant stabilizer band for separate use or in combination with surgical bra as a stabilization support of implants following augmentation.

Colors: 

Product description:

- 1** ergonomically shaped stabilizer band
- 2** velcro fastening
- 3** dense material made from lightweight rubber technology




Size	MEASUREMENTS IN CM / EU SIZE	
	Breast circumference	Under breasts
65	77-85	63-67
70	82-90	68-72
75	87-95	73-77
80	92-100	78-82
85	97-105	83-87
90	102-110	88-92



PI unique

LIPO-PI07C00X



Use: Postoperative care after plastic-aesthetic procedures: Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors: ● ●

Product description:

- 1** butterfly-shaped 5-positions adjustable shoulder straps
- 2** reinforced bottom cup area
- 3** zipper underlain with hooks
- 4** soft elastic band with logo




Size	MEASUREMENTS IN CM / EU SIZE				
	Under breasts	Breast circumference			
		Cup A	Cup B	Cup C	Cup D
65	63-67	77-79	79-81	81-83	83-85
70	68-72	82-84	84-86	86-88	88-90
75	73-77	87-89	89-91	91-93	93-95
80	78-82	92-94	94-96	96-98	98-100
85	83-87	97-99	99-101	101-103	103-105
90	88-92	102-104	104-106	106-108	108-110



PI premium

LIPO-PI20V00X



Use: Fashionable bra PI premium is specifically designed for women who have already finished their post-operative recovery period after breast implants operation.*

* The length of post-operative care is decided individually by a surgeon. PI premium bra is suitable straight after the recovery is finished.

PI premium is not a Medical device.


Colors: ● ●

Product description:

- 1** straps with an adjustable buckle
- 2** closing on the back as normal bra
- 3** elegant design with lace and golden accessories – nickel-free
- 4** heat moulded seamless cups



Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Breast circumference	
		Cup A - B	Cup C - D
65	63-67	77-81	81-85
70	68-72	82-86	86-90
75	73-77	87-91	91-95
80	78-82	92-96	96-100
85	83-87	97-101	101-105
90	88-92	102-106	106-110
95	93-97	107-111	111-115



WHEN ELEGANCE MEETS COMPRESSION

„Based on my own experience after augmentation, I know that buying an elegant bra that keeps my new breasts safe is a challenging task. PI premium has all that every woman expects from a bra. It is elegant and comfortable, and thanks to its compression, also helps to fix implants in the right place. I'm very proud of our PI premium, and I'm looking forward to your feedback.“

*Alena Hradilova,
Co-Founder of LIPOELASTIC a.s.*



EST 2002

lipoelastic®

MEDICAL PRODUCTS




LIPOELASTIC

FEMALE GARMENTS

VH Comfort

LIPO-VH00C00C



Use: Postoperative care after cesarean delivery and after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

Colors: 

Product description:

- 1** soft elastic band with logo
- 2** side zipper underlain with hooks
- 3** the abdomen area is reinforced with double layer of fabric for optimal compression
- 4** 3 positions of hooks-and-eyes closure




Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Waist	Hips
XS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128
3XL	107-115	104-112	129-137



VH

LIPO-VH00W00C



Use: Postoperative care after cesarean delivery and after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

Colors: 

Product description:

- 1** soft elastic band with logo
- 2** the abdomen area is reinforced with double layer of fabric for optimal compression
- 3** 'pull-up' design
- 4** 3 positions of hooks-and-eyes closure




Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Waist	Hips
XS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128
3XL	107-115	104-112	129-137



VH Variant

LIPO-VH00V00C



Use: Postoperative care after cesarean delivery and after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

Colors:  

Product description:

- 1** detachable straps with an adjustable buckle
- 2** 3 positions of hooks-and-eyes closure on both sides
- 3** the abdomen area is reinforced with double layer of fabric for optimal compression
- 4** 3 positions of hooks-and-eyes closure

BEST
SELLER




Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Waist	Hips
XS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128
3XL	107-115	104-112	129-137



VH low Variant

LIPO-VH00V09C



Use: Postoperative care after cesarean delivery and after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

Recommended for women with height up to 165 cm.

Colors: 

Product description:

- 1** detachable straps with an adjustable buckle
- 2** 3 positions of hooks-and-eyes closure on both sides
- 3** the abdomen area is reinforced with double layer of fabric for optimal compression




Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Waist	Hips
XS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128
3XL	107-115	104-112	129-137



VH unique Variant

LIPO-VH07V00U



Use: Postoperative care after cesarean delivery and after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

Colors:  

Product description:

- 1** detachable straps with an adjustable buckle
- 2** 3 positions of hooks-and-eyes closure on one side
- 3** the abdomen area is reinforced with double layer of fabric for optimal compression
- 4** 3 positions of hooks-and-eyes closure




Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Waist	Hips
XS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128
3XL	107-115	104-112	129-137



VH body Variant

LIPO-VH12V00C



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

Colors: 

Product description:

- 1** non-detachable straps with 4 positions of hooks-and-eyes closure
- 2** high back
- 3** 3 positions of hooks-and-eyes closure
- 4** the abdomen area is reinforced with double layer of fabric for optimal compression
- 5** 3 positions of hooks-and-eyes closure




Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Waist	Hips
XS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128
3XL	107-115	104-112	129-137



VH special Comfort

LIPO-VH01C00C




Use: Postoperative care after cesarean delivery and after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

Colors: 

Product description:

- 1** detachable straps with an adjustable buckle
- 2** both side zippers underlain with hooks on both sides for easier dressing
- 3** foam padding in the abdominal area
- 4** 3 positions of hooks-and-eyes closure

DOCTOR'S
RECOMMENDATION




Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Waist	Hips
XS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128
3XL	107-115	104-112	129-137



VH special

LIPO-VH01W00C



Use: Postoperative care after cesarean delivery and after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

Colors: 

Product description:

- 1** detachable straps with an adjustable buckle
- 2** foam padding in the abdominal area
- 3** 'pull-up' design
- 4** 3 positions of hooks-and-eyes closure




Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Waist	Hips
XS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128
3XL	107-115	104-112	129-137



VF Comfort

LIPO-VF00C00C



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, thighs, back

Colors: 

Product description:

- 1** detachable straps with an adjustable buckle
- 2** zipper underlain with hooks on both sides for easier dressing
- 3** the abdomen area is reinforced with double layer of fabric for optimal compression
- 4** the outer part of the thigh is reinforced with double layer of fabric for optimal compression
- 5** hygiene opening in the crotch for all day comfort




Size	MEASUREMENTS IN CM / EU SIZE				
	Under breasts	Waist	Hips	Thigh	Above knee
XS	61-66	58-63	84-88	48-50	32-34
S	67-72	64-69	89-94	51-55	35-39
M	73-79	70-76	95-101	56-60	40-44
L	80-88	77-85	102-110	61-65	45-49
XL	89-97	86-94	111-119	66-70	50-54
XXL	98-106	95-103	120-128	71-75	55-59
3XL	107-115	104-112	129-137	76-80	60-64



VF Variant

LIPO-VF00V00C



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, thighs, back


Colors: 

Product description:

- 1** detachable straps with an adjustable buckle
- 2** 3 positions of hooks-and-eyes closure
- 3** the abdomen area is reinforced with double layer of fabric for optimal compression
- 4** the outer part of the thigh is reinforced with double layer of fabric for optimal compression
- 5** hygiene opening in the crotch for all day comfortt




Size	MEASUREMENTS IN CM / EU SIZE				
	Under breasts	Waist	Hips	Thigh	Above knee
XS	61-66	58-63	84-88	48-50	32-34
S	67-72	64-69	89-94	51-55	35-39
M	73-79	70-76	95-101	56-60	40-44
L	80-88	77-85	102-110	61-65	45-49
XL	89-97	86-94	111-119	66-70	50-54
XXL	98-106	95-103	120-128	71-75	55-59
3XL	107-115	104-112	129-137	76-80	60-64



VF body Variant

LIPO-VF12V00C



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, thighs, back

Colors: 

Product description:


- 1** non-detachable straps with 4 positions of hooks-and-eyes closure
- 2** high back
- 3** 3 positions of hooks-and-eyes closure
- 4** the abdomen area is reinforced with double layer of fabric for optimal compression
- 5** the outer part of the thigh is reinforced with double layer of fabric for optimal compression
- 6** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE				
	Under breasts	Waist	Hips	Thigh	Above knee
XS	61-66	58-63	84-88	48-50	32-34
S	67-72	64-69	89-94	51-55	35-39
M	73-79	70-76	95-101	56-60	40-44
L	80-88	77-85	102-110	61-65	45-49
XL	89-97	86-94	111-119	66-70	50-54
XXL	98-106	95-103	120-128	71-75	55-59
3XL	107-115	104-112	129-137	76-80	60-64



VF body BBO Variant

LIPO-VF13V00C




Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Fat transfer to the buttocks area, Liposuction of the abdomen, flanks, thighs, back

Colors: 

Product description:


- 1** non-detachable straps with 4 positions of hooks-and-eyes closure
- 2** high back
- 3** 3 positions of hooks-and-eyes closure
- 4** open buttock design
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE				
	Under breasts	Waist	Hips	Thigh	Above knee
XS	61-66	58-63	84-88	48-50	32-34
S	67-72	64-69	89-94	51-55	35-39
M	73-79	70-76	95-101	56-60	40-44
L	80-88	77-85	102-110	61-65	45-49
XL	89-97	86-94	111-119	66-70	50-54
XXL	98-106	95-103	120-128	71-75	55-59
3XL	107-115	104-112	129-137	76-80	60-64



VF body BBC Variant

LIPO-VF14V00C



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Fat transfer to the buttocks area, Liposuction of the abdomen, flanks, thighs, back

Colors: 

Product description:

- 1** non-detachable straps with 4 positions of hooks-and-eyes closure
- 2** high back
- 3** 3 positions of hooks-and-eyes closure
- 4** closed buttock design (without compression on the buttocks)
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE				
	Under breasts	Waist	Hips	Thigh	Above knee
XS	61-66	58-63	84-88	48-50	32-34
S	67-72	64-69	89-94	51-55	35-39
M	73-79	70-76	95-101	56-60	40-44
L	80-88	77-85	102-110	61-65	45-49
XL	89-97	86-94	111-119	66-70	50-54
XXL	98-106	95-103	120-128	71-75	55-59
3XL	107-115	104-112	129-137	76-80	60-64



VD Comfort

LIPO-VD00C00C




Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, back

Colors: 

Product description:


- 1** detachable straps with an adjustable buckle
- 2** zipper underlain with hooks on both sides for easier dressing
- 3** the abdomen area is reinforced with double layer of fabric for optimal compression
- 4** the outer part of the thigh is reinforced with double layer of fabric for optimal compression
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
XS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44
3XL	107-115	104-112	129-137	76-80	60-64	45-47



VD Variant

LIPO-VD00V00C



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, back

Colors: 

Product description:

- 1** detachable straps with an adjustable buckle
- 2** 3 positions of hooks-and-eyes closure
- 3** the abdomen area is reinforced with double layer of fabric for optimal compression
- 4** the outer part of the thigh is reinforced with double layer of fabric for optimal compression
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
XS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44
3XL	107-115	104-112	129-137	76-80	60-64	45-47



VD body Variant

LIPO-VD12V00C




Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, back

Colors: 

Product description:


- 1** non-detachable straps with 4 positions of hooks-and-eyes closure
- 2** high back
- 3** 3 positions of hooks-and-eyes closure
- 4** the abdomen area is reinforced with double layer of fabric for optimal compression
- 5** the outer part of the thigh is reinforced with double layer of fabric for optimal compression
- 6** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
XS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44
3XL	107-115	104-112	129-137	76-80	60-64	45-47



VD body BBO Variant

LIPO-VD13V00C



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Fat transfer to the buttocks area, Liposuction of the abdomen, flanks, thighs, knees, back

Colors: 

Product description:

- 1** non-detachable straps with 4 positions of hooks-and-eyes closure
- 2** high back
- 3** 3 positions of hooks-and-eyes closure
- 4** open buttock design
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
XS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44
3XL	107-115	104-112	129-137	76-80	60-64	45-47



VD body BBC Variant

LIPO-VD14V00C



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Fat transfer to the buttocks area, Liposuction of the abdomen, flanks, thighs, knees, back

Colors: 

Product description:


- 1** non-detachable straps with 4 positions of hooks-and-eyes closure
- 2** high back
- 3** 3 positions of hooks-and-eyes closure
- 4** closed buttock design (without compression on the buttocks)
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
XS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44
3XL	107-115	104-112	129-137	76-80	60-64	45-47



VD special Comfort

LIPO-VD01C00C




Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, back

Colors: 

Product description:

- 1** detachable straps with an adjustable buckle
- 2** zipper underlain with hooks on both sides for easier dressing
- 3** the abdomen area is reinforced with double layer of fabric for optimal compression
- 4** the outer and inner part of the thigh is reinforced with double layer of fabric for optimal compression
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
XS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44
3XL	107-115	104-112	129-137	76-80	60-64	45-47



VB Comfort

LIPO-VB00C00C



Use: Disorders of the lymphatic system of the lower limbs, in patients suffering from lipedema of the lower limbs and during postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, calves, back

Colors: 

Product description:


- 1** detachable straps with an adjustable buckle
- 2** zipper underlain with hooks on both sides for easier dressing
- 3** garment finished with a clean-cut
- 4** the abdomen area is reinforced with double layer of fabric for optimal compression
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
XS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44
3XL	107-115	104-112	129-137	76-80	60-64	45-47



VB Variant

LIPO-VB00V00C



Use: Disorders of the lymphatic system of the lower limbs, in patients suffering from lipedema of the lower limbs and during postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, calves, back

Colors: 

Product description:

- 1** detachable straps with an adjustable buckle
- 2** 3 positions of hooks-and-eyes closure
- 3** the abdomen area is reinforced with double layer of fabric for optimal compression
- 4** garment finished with a clean-cut
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
XS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44
3XL	107-115	104-112	129-137	76-80	60-64	45-47



VBfLg Variant

📏 **LOW:** LIPO-VBfLg00V09C 🔗 **HIGH:** LIPO-VBfLg00V10C



Use: Lower limb swelling; disorders of the lymphatic system of the lower limbs, in patients suffering from lipedema of the lower limbs, and during postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, calves, back

LOW: Recommended for body height between 155 and 165 cm
HIGH: Recommended for body height between 165 and 175 cm

Other body sizes can be done as made to measure.


Colors: ● ●



Product description:


- 1** 3 positions of hooks-and-eyes closure
- 2** graduated compression
- 3** detachable straps with an adjustable buckle
- 4** garment finished with a clean-cut
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
XS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44
3XL	107-115	104-112	129-137	76-80	60-64	45-47



TF Comfort

LIPO-TF00C00C



Use: Postoperative care after plastic-aesthetic procedures: Liposuction of the abdomen, buttocks, hip area, flanks, thighs

Colors: 

Product description:

- 1** both side zippers underlain with hooks on both sides for easier dressing
- 2** the outer part of the thigh is reinforced with double layer of fabric for optimal compression
- 3** waistline finished with hemming technology and secured with a fine seam
- 4** hygiene opening in the crotch for all day comfort



Size	MEASUREMENTS IN CM / EU SIZE			
	Waist	Hips	Thigh	Above knee
XS	58-63	84-88	48-50	32-34
S	64-69	89-94	51-55	35-39
M	70-76	95-101	56-60	40-44
L	77-85	102-110	61-65	45-49
XL	86-94	111-119	66-70	50-54
XXL	95-103	120-128	71-75	55-59
3XL	104-112	129-137	76-80	60-64

TF Variant

LIPO-TF00V00C




Use: Postoperative care after plastic-aesthetic procedures: Liposuction of the abdomen, buttocks, hip area, flanks, thighs

Colors: 

Product description:


- 1** 3 positions of hooks-and-eyes closure
- 2** waistline finished with using hemming technology and secured with a fine seam
- 3** the outer part of the thigh is reinforced with double layer of fabric for optimal compression
- 4** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE			
	Waist	Hips	Thigh	Above knee
XS	58-63	84-88	48-50	32-34
S	64-69	89-94	51-55	35-39
M	70-76	95-101	56-60	40-44
L	77-85	102-110	61-65	45-49
XL	86-94	111-119	66-70	50-54
XXL	95-103	120-128	71-75	55-59
3XL	104-112	129-137	76-80	60-64



TF leggings

LIPO-TF15W00C



Use: For use in the symptoms of lower limbs cellulite and during postoperative care after plastic-aesthetic procedures: Liposuctions of the abdomen, buttocks, hip area, flanks, thighs

Colors: 

Product description:

- 1** soft elastic band with logo
- 2** no hygiene opening



DOCTOR'S
RECOMMENDATION



Size	MEASUREMENTS IN CM / EU SIZE			
	Waist	Hips	Thigh	Above knee
XS	58-63	84-88	48-50	32-34
S	64-69	89-94	51-55	35-39
M	70-76	95-101	56-60	40-44
L	77-85	102-110	61-65	45-49
XL	86-94	111-119	66-70	50-54
XXL	95-103	120-128	71-75	55-59
3XL	104-112	129-137	76-80	60-64

TD Comfort

LIPO-TD00C00C



Use: Postoperative care after plastic-aesthetic procedures: Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees

Colors: 

Product description:

- 1** zipper underlain with hooks on both sides for easier dressing
- 2** the outer part of the thigh is reinforced with double layer of fabric for optimal compression
- 3** waistline finished with hemming technology and secured with a fine seam
- 4** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE				
	Waist	Hips	Thigh	Above knee	Calf
XS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44
3XL	104-112	129-137	76-80	60-64	45-47



TD Variant

LIPO-TD00V00C



Use: Postoperative care after plastic-aesthetic procedures: Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees

Colors: 

Product description:

- 1** 3 positions of hooks-and-eyes closure
- 2** the outer part of the thigh is reinforced with double layer of fabric for optimal compression
- 3** waistline finished with hemming technology and secured with a fine seam
- 4** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE				
	Waist	Hips	Thigh	Above knee	Calf
XS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44
3XL	104-112	129-137	76-80	60-64	45-47



TD leggings

LIPO-TD15W00C





Use: For symptoms of lower limb cellulite; and during postoperative care after plastic-aesthetic procedures: Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees

Colors: 

Product description:


- 1** soft elastic band with logo
- 2** no hygiene opening



Size	MEASUREMENTS IN CM / EU SIZE				
	Waist	Hips	Thigh	Above knee	Calf
XS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44
3XL	104-112	129-137	76-80	60-64	45-47

TB Comfort

LIPO-TB00C00C




Use: Disorders of the lymphatic system of the lower limbs, in patients suffering from lipedema of the lower limbs and during postoperative care after plastic-aesthetic procedures: Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, calves

Colors: 

Product description:


- 1** waistline finished with hemming technology and secured with a fine seam
- 2** garment finished with a clean-cut
- 3** zipper underlain with hooks on both sides for easier dressing
- 4** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE				
	Waist	Hips	Thigh	Above knee	Calf
XS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44
3XL	104-112	129-137	76-80	60-64	45-47



TB Variant

LIPO-TB00V00C



Use: Disorders of the lymphatic system of the lower limbs, for patients suffering from lipedema of the lower limbs and during postoperative care after plastic-aesthetic procedures: Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, calves

Colors: 

Product description:


- 1** 3 positions of hooks-and-eyes closure
- 2** garment finished with a clean-cut
- 3** waistline finished with hemming technology and secured with a fine seam
- 4** hygiene opening in the crotch for all day comfort



Size	MEASUREMENTS IN CM / EU SIZE				
	Waist	Hips	Thigh	Above knee	Calf
XS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44
3XL	104-112	129-137	76-80	60-64	45-47

TB leggings

LIPO-TB15W00C



Use: The symptoms of lower limbs cellulite; lower limb swelling; disorders of the lymphatic system of the lower limbs, in patients suffering from lipedema of the lower limbs, and during postoperative care after plastic-aesthetic procedures: Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, calves

Colors: 

Product description:

- 1** soft elastic band with logo
- 2** no hygiene opening





BEST
SELLER



Size	MEASUREMENTS IN CM / EU SIZE				
	Waist	Hips	Thigh	Above knee	Calf
XS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44
3XL	104-112	129-137	76-80	60-64	45-47

TBfLg Variant

✓ **LOW:** LIPO-TBfLg00V09C ↗ **HIGH:** LIPO-TBfLg00V10C



Use: Lower limb swelling; disorders of the lymphatic system of the lower limbs, in patients suffering from lipedema of the lower limbs, and during postoperative care after plastic-aesthetic procedures: Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, calves

LOW: Recommended for body height between 155 and 165 cm

HIGH: Recommended for body height between 165 and 175 cm

Other body sizes can be done as made to measure.

Colors: ● ●



Product description:

- 1** 3 positions of hooks-and-eyes closure
- 2** graduated compression
- 3** garment finished with a clean-cut
- 4** waistline finished with hemming technology and secured with a fine seam
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE				
	Waist	Hips	Thigh	Above knee	Calf
XS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44
3XL	104-112	129-137	76-80	60-64	45-47




BEST
SELLER



TBfLg leggings

↙ **LOW:** LIPO-TBfLg28W09C ↗ **HIGH:** LIPO-TBfLg28W10C



Use: Lower limb swelling; disorders of the lymphatic system of the lower limbs, in patients suffering from lipedema of the lower limbs, and during postoperative care after plastic-aesthetic procedures: Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, calves

LOW: Recommended for body height between 155 and 165 cm

HIGH: Recommended for body height between 165 and 175 cm

Colors: 

Product description:

- 1** soft elastic band with logo
- 2** garment finished with a clean-cut
- 3** no hygiene opening



Size	MEASUREMENTS IN CM				
	Waist	Hips	Thigh	Above knee	Calf
XS	58-63	84-88	48-50	32-34	27-29
XS+	58-63	84-88	48-50	32-34	30-32
S	64-69	89-94	51-55	35-39	30-32
S+	64-69	89-94	51-55	35-39	33-35
M	70-76	95-101	56-60	40-44	33-35
M+	70-76	95-101	56-60	40-44	36-38
L	77-85	102-110	61-65	45-49	36-38
L+	77-85	102-110	61-65	45-49	39-41
XL	86-94	111-119	66-70	50-54	39-41
XL+	86-94	111-119	66-70	50-54	42-44
XXL	95-103	120-128	71-75	55-59	42-44
XXL+	95-103	120-128	71-75	55-59	45-47
3XL	104-112	129-137	76-80	60-64	45-47
3XL+	104-112	129-137	76-80	60-64	48-50



NEW



ACTIVE leggings

LIPO-ACTIVE



Use: The symptoms of lower limbs cellulite; lower limbs swelling; disorders of the lymphatic system of the lower limbs, in patients suffering from lipedema of the lower limbs, as well as during postoperative care after cesarean delivery and in postoperative care after plastic-aesthetic procedures: Liposuction of the abdomen, buttocks, hip area, flanks, thighs, knees, calves



ACTIVE leggings are suitable for everyday wear, travel, and sports.

Colors: ● ● ●

Product description:

- 1** fashionable soft elastic band
- 2** garment finished with a clean-cut
- 3** no hygiene opening

Size	MEASUREMENTS IN CM				
	Waist	Hips	Thigh	Above knee	Calf
XS	58-63	84-88	48-50	32-34	27-29
XS+	58-63	84-88	48-50	32-34	30-32
S	64-69	89-94	51-55	35-39	30-32
S+	64-69	89-94	51-55	35-39	33-35
M	70-76	95-101	56-60	40-44	33-35
M+	70-76	95-101	56-60	40-44	36-38
L	77-85	102-110	61-65	45-49	36-38
L+	77-85	102-110	61-65	45-49	39-41
XL	86-94	111-119	66-70	50-54	39-41
XL+	86-94	111-119	66-70	50-54	42-44
XXL	95-103	120-128	71-75	55-59	42-44
XXL+	95-103	120-128	71-75	55-59	45-47
3XL	104-112	129-137	76-80	60-64	45-47
3XL+	104-112	129-137	76-80	60-64	48-50







EST 2002

lipoelastic[®]

MEDICAL PRODUCTS

LIPOELASTIC

ARM SLEEVES

AP Variant

LIPO-AP00V00C



Use: Postoperative care after plastic-aesthetic procedures: Brachioplasty, Liposuction of the arms, armpits, back

Colors: 

Product description:

- 1** 3 positions of hooks-and-eyes closure
- 2** offers special padding foam in the armpit area to prevent discomfort
- 3** sleeve finished with a clean-cut



Size	MEASUREMENTS IN CM / EU SIZE	
	Arm circumference	Under elbow circumference
XS	26-27	21-22
S	28-29	23-24
M	30-31	25-26
L	32-34	27-28
XL	35-37	29-31
XXL	38-40	32-34
3XL	41-43	35-37



AP long Variant

LIPO-AP00V07C



Use: Disorders of the lymphatic system of the upper limbs, in patients suffering from lipedema of the upper limbs and during postoperative care after plastic-aesthetic procedures: Brachioplasty, Liposuction of the arms, armpits, back

Colors: 

Product description:

- 1** 3 positions of hooks-and-eyes closure
- 2** offers special padding foam in the armpit area to prevent discomfort
- 3** long sleeve finished with a clean-cut

BEST
SELLER



Size	MEASUREMENTS IN CM / EU SIZE	
	Arm circumference	Under elbow circumference
XS	26-27	21-22
S	28-29	23-24
M	30-31	25-26
L	32-34	27-28
XL	35-37	29-31
XXL	38-40	32-34
3XL	41-43	35-37



APf Variant

LIPO-APf00V07C



Use: For female patients for use in lymphedema and lipedema of the upper limbs and during postoperative care after plastic-aesthetic procedures: Brachioplasty, Liposuction of the arms, back

Colors:  

Product description:

- 1** 3 positions of hooks-and-eyes closure
- 2** fixation on finger is finished with a clean-cut
- 3** breast area finished with a clean-cut
- 4** seamless band finished with hemming technology and secured with a fine seam




Size	MEASUREMENTS IN CM / EU SIZE		
	Arm circumference	Under elbow circumference	Under breasts
XS	26-27	21-22	61-66
S	28-29	23-24	67-72
M	30-31	25-26	73-79
L	32-34	27-28	80-88
XL	35-37	29-31	89-97
XXL	38-40	32-34	98-106
3XL	41-43	35-37	107-115



AS Variant

LIPO-AS00V00C



Use: Postoperative care after plastic-aesthetic procedures: Brachioplasty, Liposuction of the arms

Colors: 

Product description:

- 1** 3 positions of hooks-and-eyes back closure
- 2** sleeve finished with a clean-cut




Size	MEASUREMENTS IN CM / EU SIZE	
	Arm circumference	Under elbow circumference
XS	26-27	21-22
S	28-29	23-24
M	30-31	25-26
L	32-34	27-28
XL	35-37	29-31
XXL	38-40	32-34
3XL	41-43	35-37



AS long Variant

LIPO-AS00V07C



Use: Disorders of the lymphatic system of the upper limbs, in patients suffering from lipedema of the upper limbs and during postoperative care after plastic-aesthetic procedures: Brachioplasty, Liposuction of the arms

Colors: 

Product description:

- 1** 3 positions of hooks and eyes back closure
- 2** long sleeve finished with a clean-cut



Size	MEASUREMENTS IN CM / EU SIZE	
	Arm circumference	Under elbow circumference
XS	26-27	21-22
S	28-29	23-24
M	30-31	25-26
L	32-34	27-28
XL	35-37	29-31
XXL	38-40	32-34
3XL	41-43	35-37



EST 2002

lipoelastic®

MEDICAL PRODUCTS



LIPOELASTIC

FEMALE VESTS
AND BODY SUITS

MT smooth Variant

LIPO-MT08V00R




Use: Disorders of the lymphatic system of the upper limbs and for use during postoperative care after plastic-aesthetic procedures: Brachioplasty, Liposuction of the arms, back, Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors: 

Product description:

- 1** soft elastic band with logo
- 2** moulded cups
- 3** 3 positions of hooks-and-eyes closure
- 4** sleeve finished with a clean-cut




Size	MEASUREMENTS IN CM / EU SIZE			
	Breast circumference	Under breasts	Waist	Arm
XS	72-79	61-66	58-63	26-27
S	80-87	67-72	64-69	28-29
M	88-95	73-79	70-76	30-31
L	96-103	80-88	77-85	32-34
XL	104-111	89-97	86-94	35-37
XXL	112-120	98-106	95-103	38-40
3XL	121-129	107-115	104-112	41-43



MTf smooth Comfort

LIPO-MTF08C00C



Use: Disorders of the lymphatic system of the upper limbs, in patients suffering from lipedema of the upper limbs and during postoperative care after plastic-aesthetic procedures: Brachioplasty, Liposuction of the arms, back

Colors: 

Product description:

- 1** front zipper underlain with hooks
- 2** soft elastic band with logo
- 3** sleeve fixation on finger
- 4** fixation on finger is finished with a clean-cut




BEST
SELLER

Size	MEASUREMENTS IN CM / EU SIZE			
	Breast circumference	Under breasts	Waist	Arm
XS	72-79	61-66	58-63	26-27
S	80-87	67-72	64-69	28-29
M	88-95	73-79	70-76	30-31
L	96-103	80-88	77-85	32-34
XL	104-111	89-97	86-94	35-37
XXL	112-120	98-106	95-103	38-40
3XL	121-129	107-115	104-112	41-43



MH Comfort

LIPO-MH00C00C



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Brachioplasty, Liposuction of the arms, abdomen, flanks, back

Colors: 

Product description:

- 1** front zipper underlain with hooks
- 2** the abdomen area is reinforced with double layer of fabric for optimal compression
- 3** sleeve finished with a clean-cut
- 4** 3 positions of hooks-and-eyes closure




Size	MEASUREMENTS IN CM / EU SIZE				
	Breast circumference	Under breasts	Waist	Hips	Arm
XS	72-79	61-66	58-63	84-88	26-27
S	80-87	67-72	64-69	89-94	28-29
M	88-95	73-79	70-76	95-101	30-31
L	96-103	80-88	77-85	102-110	32-34
XL	104-111	89-97	86-94	111-119	35-37
XXL	112-120	98-106	95-103	120-128	38-40
3XL	121-129	107-115	104-112	129-137	41-43



MH special Comfort

LIPO-MH01C00R



Use: Postoperative care after cesarean delivery and after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back, Mammoplasty (Augmentation, Breast reduction, Breast reconstruction), Mastopexy

Colors: 

Product description:

- 1** straps with 4 positions of hooks-and-eyes closure
- 2** moulded cups
- 3** front zipper underlain with hooks
- 4** 3 positions of hooks-and-eyes closure




Size	MEASUREMENTS IN CM / EU SIZE			
	Breast circumference	Under breasts	Waist	Hips
XS	72-79	61-66	58-63	84-88
S	80-87	67-72	64-69	89-94
M	88-95	73-79	70-76	95-101
L	96-103	80-88	77-85	102-110
XL	104-111	89-97	86-94	111-119
XXL	112-120	98-106	95-103	120-128
3XL	121-129	107-115	104-112	129-137



MHF Comfort

LIPO-MHF00C00C




Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Brachioplasty, Liposuction of the arms, abdomen, buttocks, hip area, flanks, thighs, back

Colors: 

Product description:

- 1** front zipper underlain with hooks
- 2** the abdomen area is reinforced with double layer of fabric for optimal compression
- 3** hygiene opening in the crotch for all day comfort



Size	MEASUREMENTS IN CM / EU SIZE						
	Breast circumference	Under breasts	Waist	Hips	Thigh	Above knee	Arm
XS	72-79	61-66	58-63	84-88	48-50	32-34	26-27
S	80-87	67-72	64-69	89-94	51-55	35-39	28-29
M	88-95	73-79	70-76	95-101	56-60	40-44	30-31
L	96-103	80-88	77-85	102-110	61-65	45-49	32-34
XL	104-111	89-97	86-94	111-119	66-70	50-54	35-37
XXL	112-120	98-106	95-103	120-128	71-75	55-59	38-40
3XL	121-129	107-115	104-112	129-137	76-80	60-64	41-43

MHD Comfort

LIPO-MHD00C00C




Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Brachioplasty, Liposuction of the arms, abdomen, buttocks, hip area, flanks, thighs, knees, back

Colors: 

Product description:


- 1** front zipper underlain with hooks
- 2** the abdomen area is reinforced with double layer of fabric for optimal compression
- 3** hygiene opening in the crotch for all day comfort



Size	MEASUREMENTS IN CM / EU SIZE							
	Breast circumference	Under breasts	Waist	Hips	Thigh	Above knee	Calf	Arm
XS	72-79	61-66	58-63	84-88	48-50	32-34	27-29	26-27
S	80-87	67-72	64-69	89-94	51-55	35-39	30-32	28-29
M	88-95	73-79	70-76	95-101	56-60	40-44	33-35	30-31
L	96-103	80-88	77-85	102-110	61-65	45-49	36-38	32-34
XL	104-111	89-97	86-94	111-119	66-70	50-54	39-41	35-37
XXL	112-120	98-106	95-103	120-128	71-75	55-59	42-44	38-40
3XL	121-129	107-115	104-112	129-137	76-80	60-64	45-47	41-43

MGF Variant

LIPO-MGF00V00R



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, back, buttocks, hip area, flanks, thighs, Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors:

Product description:

- 1** straps with 4 positions of hooks-and-eyes closure
- 2** moulded cups
- 3** garment finished with a clean-cut
- 4** 3 positions of hooks-and-eyes closure
- 5** hygiene opening in the crotch for all day comfort



Size	MEASUREMENTS IN CM / EU SIZE					
	Breast circumference	Under breasts	Waist	Hips	Thigh	Above knee
XS	72-79	61-66	58-63	84-88	48-50	32-34
S	80-87	67-72	64-69	89-94	51-55	35-39
M	88-95	73-79	70-76	95-101	56-60	40-44
L	96-103	80-88	77-85	102-110	61-65	45-49
XL	104-111	89-97	86-94	111-119	66-70	50-54
XXL	112-120	98-106	95-103	120-128	71-75	55-59
3XL	121-129	107-115	104-112	129-137	76-80	60-64

MGD Variant

LIPO-MGD00V00R




Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, back, buttocks, hip area, flanks, thighs, knees, Mammoplasty (Augmentation, Breast Reduction, Breast Reconstruction), Mastopexy

Colors: 

Product description:

- 1** straps with 4 positions of hooks-and-eyes closure
- 2** moulded cups
- 3** 3 positions of hooks-and-eyes closure
- 4** garment finished with a clean-cut
- 5** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE						
	Breast circumference	Under breasts	Waist	Hips	Thigh	Above knee	Calf
XS	72-79	61-66	58-63	84-88	48-50	32-34	27-29
S	80-87	67-72	64-69	89-94	51-55	35-39	30-32
M	88-95	73-79	70-76	95-101	56-60	40-44	33-35
L	96-103	80-88	77-85	102-110	61-65	45-49	36-38
XL	104-111	89-97	86-94	111-119	66-70	50-54	39-41
XXL	112-120	98-106	95-103	120-128	71-75	55-59	42-44
3XL	121-129	107-115	104-112	129-137	76-80	60-64	45-47







EST 2002

lipoelastic®


MEDICAL PRODUCTS

LIPOELASTIC

MALE GARMENTS

VHmS Comfort

LIPO-VHmS00C00R



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

Colors: 

Product description:

- 1** soft elastic band with logo
- 2** front zipper underlain with hooks
- 3** hygiene opening in the crotch for all day comfort




Size	MEASUREMENTS IN CM / EU SIZE		
	Under chest circumference	Waist	Hips
XS	82-87	75-79	82-85
S	88-92	80-86	86-90
M	93-98	87-93	91-97
L	99-105	94-101	98-104
XL	106-113	102-110	105-112
XXL	114-121	111-119	113-120
3XL	122-130	120-128	121-128



VHmm Comfort

LIPO-VHmm00C00R



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, back

Colors: 

Product description:

- 1** soft elastic band with logo
- 2** zipper underlain with hooks
- 3** garment finished with a clean-cut
- 4** hygiene opening in the crotch for all day comfort




Size	MEASUREMENTS IN CM / EU SIZE		
	Under chest circumference	Waist	Hips
XS	82-87	75-79	82-85
S	88-92	80-86	86-90
M	93-98	87-93	91-97
L	99-105	94-101	98-104
XL	106-113	102-110	105-112
XXL	114-121	111-119	113-120
3XL	122-130	120-128	121-128



VFm Comfort

LIPO-VFm00C00R



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, buttocks, hip area, flanks, thighs, back

Colors: 

Product description:

- 1** detachable straps with an adjustable buckle
- 2** zipper underlain with hooks on both sides for easier dressing
- 3** garment finished with a clean-cut
- 4** hygiene opening in the crotch for all day comfort




Size	MEASUREMENTS IN CM / EU SIZE				
	Under chest circumference	Waist	Hips	Thigh	Above knee
XS	82-87	75-79	82-85	47-50	31-34
S	88-92	80-86	86-90	51-55	35-39
M	93-98	87-93	91-97	56-60	40-44
L	99-105	94-101	98-104	61-65	45-49
XL	106-113	102-110	105-112	66-70	50-54
XXL	114-121	111-119	113-120	71-75	55-59
3XL	122-130	120-128	121-128	76-80	60-64



TDm Variant

LIPO-TDm00V00R



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the lower abdomen, lower back, hips, buttocks, thighs, knees

Colors: 

Product description:


- 1** 3 positions of hooks-and-eyes closure
- 2** garment finished with a clean-cut
- 3** hygiene opening in the crotch for all day comfort

Size	MEASUREMENTS IN CM / EU SIZE				
	Waist	Hips	Thigh	Above knee	Calf
XS	75-79	82-85	47-50	31-34	30-32
S	80-86	86-90	51-55	35-39	33-35
M	87-93	91-97	56-60	40-44	36-38
L	94-101	98-104	61-65	45-49	39-41
XL	102-110	105-112	66-70	50-54	42-44
XXL	111-119	113-120	71-75	55-59	45-47
3XL	120-128	121-128	76-80	60-64	48-50



MTmS Comfort

LIPO-MTmS00C00R



Use: Postoperative care after chest and gynecomastia surgery and after plastic-aesthetic procedures: Liposuction of the back

Colors:  

Product description:

- 1** straps with 3 positions of hooks-and-eyes closure
- 2** front zipper underlain with hooks
- 3** soft elastic band with logo




Size	MEASUREMENTS IN CM / EU SIZE	
	Chest circumference	Waist
XXS	80-86	66-74
XS	87-91	75-79
S	92-96	80-86
M	97-102	87-93
L	103-109	94-101
XL	110-117	102-110
XXL	118-125	111-119
3XL	126-133	120-128



MTm Comfort

LIPO-MTm00C00R



Use: Postoperative care after chest and gynecomastia surgery and after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, back

Colors: 

Product description:


- 1** straps with 3 positions of hooks-and-eyes closure
- 2** front zipper underlain with hooks
- 3** soft elastic band with logo

Size	MEASUREMENTS IN CM / EU SIZE	
	Chest circumference	Waist
XS	87-91	75-79
S	92-96	80-86
M	97-102	87-93
L	103-109	94-101
XL	110-117	102-110
XXL	118-125	111-119
3XL	126-133	120-128



MTmL Comfort

LIPO-MTmL00C00R



Use: Postoperative care after chest and gynecomastia surgery and after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

Colors: 

Product description:

- 1** front zipper underlain with hooks
- 2** straps with 3 positions of hooks-and-eyes closure




Size	MEASUREMENTS IN CM / EU SIZE	
	Chest circumference	Waist
XS	87-91	75-79
S	92-96	80-86
M	97-102	87-93
L	103-109	94-101
XL	110-117	102-110
XXL	118-125	111-119
3XL	126-133	120-128



MTmL Variant

LIPO-MTmL00V00R



Use: Postoperative care after chest and gynecomastia surgery and after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

Colors: 

Product description:

- 1** 3 positions of hooks-and-eyes closure
- 2** straps with 3 positions of hooks-and-eyes closure




Size	MEASUREMENTS IN CM / EU SIZE	
	Chest circumference	Waist
XS	87-91	75-79
S	92-96	80-86
M	97-102	87-93
L	103-109	94-101
XL	110-117	102-110
XXL	118-125	111-119
3XL	126-133	120-128



MGmm Comfort

LIPO-MGmm00C00R



Use: Postoperative care after chest surgery
- gynecomastia and after plastic-aesthetic
procedures: Abdominoplasty, Liposuction of the
abdomen, back, flanks, buttocks, hip area, thighs

Colors: 

Product description:

- 1** straps with 3 positions of hooks-and-eyes closure
- 2** front zipper underlain with hooks
- 3** garment finished with a clean-cut
- 4** hygiene opening in the crotch for all day comfort




Size	MEASUREMENTS IN CM / EU SIZE			
	Chest circumference	Waist	Hips	Thigh
XS	87-91	75-79	82-85	47-50
S	92-96	80-86	86-90	51-55
M	97-102	87-93	91-97	56-60
L	103-109	94-101	98-104	61-65
XL	110-117	102-110	105-112	66-70
XXL	118-125	111-119	113-120	71-75
3XL	126-133	120-128	121-128	76-80



MGm long Variant

LIPO-MGm00V08R



Use: Postoperative care after plastic-aesthetic procedures: Brachioplasty, Abdominoplasty, Liposuction of the arms, abdomen, back, buttocks, hip area, flanks, thighs

Colors: 

Product description:

- 1** 3 positions of hooks-and-eyes closure
- 2** garment finished with a clean-cut
- 3** hygiene opening in the crotch for all day comfort



Size	MEASUREMENTS IN CM / EU SIZE			
	Chest circumference	Waist	Hips	Thigh
XS	87-91	75-79	82-85	47-50
S	92-96	80-86	86-90	51-55
M	97-102	87-93	91-97	56-60
L	103-109	94-101	98-104	61-65
XL	110-117	102-110	105-112	66-70
XXL	118-125	111-119	113-120	71-75
3XL	126-133	120-128	121-128	76-80



LIPOELASTIC,
TRUST THE EXPERTS
IN COMPRESSION
THERAPY.




LIPOELASTIC

ABDOMINAL
BINDERS AND
ACCESSORIES

KPress

24 cm: LIPO-KP16F03X

32 cm: LIPO-KP16F05X



Use: Postoperative care after cesarean delivery and for use by patients regardless of gender for use following postoperative care following plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

TIP:

Abdominal binders can be worn independently. However, for greater comfort and optimal results, they can also be combined with a compression post-operative garment by wearing it on the outside.

Colors: 

Product description:

- 1** flexible strap
- 2** separate velcro panels are designed for easy adjustment and closure
- 3** dense, tight material

UNISEX



KPress
Binder height
24 cm


Size	MEASUREMENTS IN CM / EU SIZE	
	Waist	Binder height
S-M	63 – 76	24, 32
L-XL	77 – 94	24, 32
XXL-3XL	95 – 112	24, 32
4XL-5XL	113 – 130	24, 32



KPlus

24 cm: LIPO-KP19F03X

32 cm: LIPO-KP19F05X



Use: Postoperative care after cesarean delivery and for use by patients regardless of gender for use in postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

TIP: Abdominal binders can be worn independently. However, for greater comfort and optimal results, they can also be combined with a compression post-operative garment by wearing it on the outside.

Colors:   

Product description:

- 1** flexible strap
- 2** separate velcro panels are designed for easy adjustment and closure
- 3** dense, tight material



KPlus
Binder height
24 cm 3 panels
32 cm 4 panels


Size	MEASUREMENTS IN CM / EU SIZE	
	Waist	Binder height
XXS-XS	53 - 62	24, 32
S-M	63 - 76	24, 32
L-XL	77 - 94	24, 32
XXL-3XL	95 - 112	24, 32
4XL-5XL	113 - 130	24, 32



KP soft

28 cm: LIPO-KP03F07X

NEW



Use: Postoperative care after cesarean delivery and for use by patients regardless of gender for use in postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

TIP:

Abdominal binders can be worn independently. However, for greater comfort and optimal results, they can also be combined with a compression post-operative garment by wearing it on the outside.

Colors: 



COTTON

Product description:

- 1** separate velcro panels are designed for easy adjustment and closure
- 2** pocket for easy fastening
- 3** flexible strap
- 4** comfortable cotton material



KP soft
Binder height
28 cm




Size	MEASUREMENTS IN CM / EU SIZE	
	Waist	Binder height
S-M	63-76	28
L-XL	77-94	28
XXL-3XL	95-112	28
4XL-5XL	113-132	28

KP special

23 cm: LIPO-KP01F02X

31 cm: LIPO-KP01F04X



Use: Postoperative care after cesarean delivery and for use by patients regardless of gender for use in postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, flanks, back

TIP: Abdominal binders can be worn independently. However, for greater comfort and optimal results, they can also be combined with a compression post-operative garment by wearing it on the outside.

Colors: ○

Product description:

- 1** flexible strap
- 2** separate velcro panels are designed for easy adjustment and closure
- 3** breathable and soft material

Size	MEASUREMENTS IN CM / EU SIZE	
	Waist	Binder height
XS	58-63	23, 31
S	64-69	23, 31
M	70-76	23, 31
L	77-85	23, 31
XL	86-94	23, 31
XXL	95-103	23, 31
3XL	104-112	23, 31

UNISEX




KP special
Binder height
23 cm



KPG

LIPO-KPG00F00X



Use: Postoperative care after gynecomastia surgery in men; in women after liposuction of the abdomen

Colors: 

Product description:

- 1** flexible strap
- 2** velcro fastening
- 3** breathable and soft material


UNISEX



Size	MEASUREMENTS IN CM / EU SIZE
	Circumference chest*/waist**
XS-S	72 - 87
M-L	88 - 103
XL-XXL	104 - 120


*when using on chest, measure the chest circumference

**when using on abdominal part, measure the waist circumference



KPad

LIPO-KPaD00W00X



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen

Colors:  



COTTON

Product description:


Unisex foam insertion that can be worn underneath abdominal binders or any compression garment covering abdominal area.



Size	MEASUREMENTS IN CM / EU SIZE
	Waist
S-M	63 – 76
L-XL	77 – 94
XXL-3XL	95 – 112

KPad 01

LIPO-KPaD01W00X



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen, hips

Colors:  



COTTON


Product description:

Unisex foam insertion that can be worn underneath abdominal binders or any compression garment covering abdomen, hips and thighs.



KPad 02

LIPO-KPaD02W00X



Use: Postoperative care after plastic-aesthetic procedures: Liposuction of the thighs

Colors:  



COTTON


Product description:

Unisex foam insertion that can be worn underneath abdominal binders or any compression garment covering abdomen and hips.



KPad 03

LIPO-KPaD03W00X



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the hip area

Colors:  



COTTON


Product description:

Unisex foam insertion that can be worn underneath abdominal binders or any compression garment covering abdomen, hips and thighs.



KPad 04

LIPO0096-04



Use: Postoperative care after plastic-aesthetic procedures: Abdominoplasty, Liposuction of the abdomen

Colors: 

Product description:


Unisex abdominal board insertion that should be worn underneath the abdominal binder or any compression garment covering the operated area.

This product is not a Medical device.



KPad 4D

LIPO0096-4D



Use: Postoperative care after plastic-aesthetic procedures: 360 Abdominoplasty, High Definition Liposuction, Liposuction of the abdomen

Colors: 

Product description:

A removable firm pad prevents bending in the abdomen area. The back part can be adjusted with Velcro fastening



Size	MEASUREMENTS IN CM / EU SIZE		
	Under breasts	Waist	Hips
XS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119

EST 2002

lipoelastic®

MEDICAL PRODUCTS




LIPOELASTIC

FACIAL GARMENTS

FM 02

LIPO-FM12F00X

NEW



Use: Postoperative care after plastic-aesthetic procedures: facial, chin, double chin, jaw, Liposuction of the chin and face, Facelift, Facial implants Insertion

Colors: 

Product description:

- 1** clean cut finishing
- 2** velcro fastening
- 3** adjustable fixation strap




Size	MEASUREMENTS IN CM / EU SIZE	
	Lower neck circumference	
XS	31-33	
S	33-35	
M	35-37	
L	37-39	
XL	39-41	



FM 03

LIPO-FM03F00X

NEW



Use: Facial and neck swelling, disorders of the lymphatic system of the face and neck and during postoperative care after plastic-aesthetic procedures: facial, chin, double chin, neck, Liposuction of the chin, double chin, face, Facial implants insertion, Otoplasty

Colors: 

Product description:

- 1** clean cut finishing
- 2** velcro fastening




Size	MEASUREMENTS IN CM / EU SIZE	
	Lower neck circumference	
XS	31-33	
S	33-35	
M	35-37	
L	37-39	
XL	39-41	



FM special

LIPO-FM01F00X



Use: Postoperative care after plastic-aesthetic procedures: facial, chin, double chin, neck; Liposuction of the chin, Chin implant insertion, Otoplasty

Colors: 

Product description:

1 velcro fastening




Size	MEASUREMENTS IN CM / EU SIZE	
	Head circumference (chin - top of the head)	
XS	50-53	
S	54-58	
M	59-63	
L	64-68	
XL	69-73	



FM extra

LIPO-FM02F00X



Use: Postoperative care after plastic-aesthetic procedures: facial, chin, double chin, neck;
Liposuction of the chin, Chin implant insertion

Colors: ○

Product description:

- 1** velcro fastening
- 2** breathable material




Size	MEASUREMENTS IN CM / EU SIZE	
	Head circumference (chin - top of the head)	
XS	50-53	
S	54-58	
M	59-63	
L	64-68	
XL	69-73	



FM 01

LIPO-FM24F00X



Use: Postoperative care after plastic-aesthetic procedures: facial, double chin; Otoplasty

Colors: ○


Product description:

1 velcro fastening



PU 03

LIPO-PU03F00C



Use: Postoperative care
after otoplasty

Colors



Product description:

1 velcro fastening



Size	MEASUREMENTS IN CM / EU SIZE	
	Head circumference (around the ears)	
XS	50-52	
S	53-55	
M	56-58	
L	59-61	
XL	62-64	
XXL	65-67	



EST 2002

lipoelastic®

MEDICAL PRODUCTS




LIPOELASTIC

OTHER GARMENTS

LIPO slide


LIPO-ESA94W00X



Use: The LIPO slide is used to simply pull the postoperative compression garment onto the patient.

LIPO slide is not a Medical device

Colors: ●



LIPO extend T and V

T: LIPO-EX95V00X

V: LIPO-EX92V00X

Use: If you need to extend the standard bandage size, use LIPO extend

Garments can be extended in the abdomen area by using an extension. For the LOW version, cut the extender above the stitch line.

There are three sizes available (according to garment's size):

XS-S

M-L

XL-3XL

LIPO extend T and V is not a Medical device

Colors: ● ●




LIPO extend AP

LIPO-EX30V00X

Use: For extend fastening on AP Variant, AP long Variant and APf.

LIPO extend AP is not a Medical device

Colors: ● ●



EST 2002

lipoelastic®

MEDICAL PRODUCTS




LIPOELASTIC

ANTI-EMBOLISM
STOCKINGS

LIPOTHROMBO AD

LIPO-LTAD0700X



Use: During pre-, intra- and postoperative care of surgical procedures, in patients with swelling of the lower limbs and in long-term patients to support hemodynamics

Package 1 pair / 12 pairs

Colors:

Product description:

- 1** marked heel area
- 2** marked size
- 3** opening facilitates easy checks for proper blood circulation




Size	MEASUREMENTS IN CM / EU SIZE									
	XS	S	M	L	XL	XXL	3XL	4XL	5XL	6XL
b	19-21	21-23	23-25	25-27	27-29	29-31	30-33	34-38	39-43	44-49
d	28-35	31-38	33-41	35-44	37-47	39-49	41-51	52-62	63-75	76-90

b – measure the ankle circumference in the narrowest place
d – measure the calf circumference in the place where the edge is supposed to be



LIPOTHROMBO AG

LIPO-LTAG0700X



Use: During pre-, intra- and postoperative care of surgical procedures, in patients with swelling of the lower limbs and in long-term patients to support hemodynamics

Package 1 pair / 6 pairs

Colors: ○

Product description:

- 1** silicone rim with a ring to determine the correct size
- 2** marked size
- 3** marked heel area
- 4** opening facilitates easy checks for proper blood circulation



Size	MEASUREMENTS IN CM / EU SIZE									
	XS	S	M	L	XL	XXL	3XL	4XL	5XL	6XL
b	19-21	21-23	23-25	25-27	27-29	29-31	30-33	34-38	39-43	44-49
g	44-51	48-55	52-59	55-67	59-67	63-70	67-73	74-85	86-95	95-105

b – measure the ankle circumference in the narrowest place
g – thigh girth is measured in the place of the intended rim edge



EST 2002

lipoelastic®

MEDICAL PRODUCTS




LIPOELASTIC

SCAR TREATMENT

LIPOELASTIC sheet LOLLIPOP

F-02555-56SLE



Use: The lollipop-shaped sheet is suitable for scar treatment in the breast area, abdominal area.

Package: 1 pair

Size: 7 x 16 cm

Minimal order: 5 pair




Product description:

- maintains the natural moisture of the skin
- helps regenerate the skin
- helps merge the scar with the surrounding tissue
- prevents the keloid and hypertrophic scars
- effective to old and new scar tissue

LIPOELASTIC sheet ANCHOR

F-02555-55SLE



Use: The anchor-shaped sheet is suitable for the treatment of scars in the chest and breasts areas.

Package: 1 pair

Size: 9,8 x 29,5 cm

Minimal order: 5 pair



Product description:

- maintains the natural moisture of the skin
- helps regenerate the skin
- helps merge the scar with the surrounding tissue
- prevents the keloid and hypertrophic scars
- effective to old and new scar tissue

LIPOELASTIC sheet STRIP01

F-02555-06SLE



Use: A strip of silicone sheet can be adjusted to the size and shape of the scar by cutting.

Package: 1 piece

Size: 5 x 20 cm

Minimal order: 10 pieces

Product description:

- maintains the natural moisture of the skin
- helps regenerate the skin
- helps merge the scar with the surrounding tissue
- prevents the keloid and hypertrophic scars
- effective to old and new scar tissue



LIPOELASTIC sheet STRIP02

F-02555-08CLE


Use: A strip of silicone sheet can be adjusted to the size and shape of the scar by cutting.

Package: 1 piece

Size: 20 x 30 cm


Product description:

- maintains the natural moisture of the skin
- helps regenerate the skin
- helps merge the scar with the surrounding tissue
- prevents the keloid and hypertrophic scars
- effective to old and new scar tissue



LIPOELASTIC gel

LIPOELASTIC gel



Use: LIPOELASTIC gel is a silicone gel enriched with vitamin E for the treatment and prevention of abnormal scars (hypertrophic and keloid).

LIPOELASTIC gel can be applied to any part of the skin and is particularly suitable for treating scars on the face, neck, at joints and in skin folds.

Weight: 10 g

Ingredients:

POLYSILOXANES, VITAMIN E ± 10%

Product description:

- silicone-based gel shows unique and proven scar treatment method
- gel is effective for new and old scar tissue
- fast drying silicone gel
- enriched with vitamin E, helps to regenerate the skin

Before



After





LIPOELASTIC a.s.
Prague, Czech Republic
info@lipoelastic.com



WWW.LIPOELASTIC.COM